

SUPERCRITICAL EXTRACTION PILOT PLANT FOR SOLID SAMPLES

Automated and computerized laboratory-pilot plant for extraction of solid samples by means of supercritical

ular-type built for easy return to factory for service. Each module (feed section, Extractor module and two Separator modules) system and electronic control system. All these systems are supervised for distributed control system PC based. The security ted in each one of these modules.

Feed system

Roy CO₂ pump, 4.7 l/h, 380 bar, SS-316 hydraulic membrane, refrigerated head. Inverter for computerized control. Check- and other components.
(-10°C) for cooling CO₂ line feed and CO₂ pump head.
Roy Co-solvent pump 0.3 l/h, 100 bar, SS-316 piston head. Inverter for computerized control. system and bypass for flow measures at high pressure in Co-solvent pump.

Column System

- Furnace for CO₂ preheater. Two control actions: heating by electrical power and cooling by furnace opening and closing.
- 350 cc Head Line vessel, 400 bar, easy closure system, for solid sample. Porous plate 20 microns. Quick connectors for agreeable work.
- Furnace for control temperature of extraction operation, internal thermocouple. Two control actions: heating by electrical power and cooling by furnace opening and closing.
- Bypass system (two-three way valves) for cleaning procedures.
- Pneumatic security valve put into operation by pressure control system.
- Pressure control system based on micrometric regulation servocontrolled valve. High precision in pressure control and fast response. Maximum pressure 340 bar.
- Rupture discs, check valves, filters and other components.

Three Separators

- 40 cc Head Line vessel, 400 bar, easy closure system, for extracts collection. Valve for sample.
- Furnace for control temperature of separation operation, internal thermocouple. Two control actions: heating by electrical power and cooling by furnace opening and closing.
- Pressure control system based on micrometric regulation servocontrolled valve. High precision in pressure control and fast response. Max pressure 220 bar in Separator 1, 120 bar in Separator 2.
- Separator 3 at atmospheric pressure and MFM for CO₂ flow measurement.

Distributed control system

ontrol systems are linked with PC computer by means of Process@ ote control with digital communications. The system can be controlled matically.

are allows the operator to design automatic procedures for the al independent safety levels: automatic switch-off in case of any re and temperature security systems, all that based on electronic or ces and independent of PC.

Test

be tested during 24 hours at 360 bar closing (except MFM and be tested during 4 hours at 340 bar in Extrac-

