

SPECIFICATIONS

HPVA II

PHYSICAL

Height	88.9 cm (35 in.)
Width	50.8 cm (20 in.)
Depth	50.8 cm (20 in.)
Weight	27.2 kg (60 lbs.)



ELECTRICAL

Voltage	100 - 240 VAC.
Frequency	50 - 60 Hz.

ENVIRONMENT

Temperature	10 to 45 °C (50 to 113 °F.) 10 to 55 °C (14 to 131 °F.)
-------------	--

COMPUTER REQUIREMENTS

Processor	Pentium IV 2.8 GHz processor.
Memory	512 MB RAM.
Hard Dish Space	80 gigabytes.
Monitor	SVGA.
Media Drive	Writable CD or DVD.
Interface	Windows XP Professional.
Software	Windows 7. Microsoft Excel - 2002 or higher.

Due to continuous improvements, specifications are subject to change without notice.