

Specifications

Storage Capacity

47 liters

Recommended Cryogens

Liquid nitrogen

Liquid argon

Electrical

Voltage: 100/110 220/240 VAC

Frequency: 50 to 60 Hz

Power: 100 VA

Weights

Storage Dewar (empty): 16.3 kg (36 lbs)

Control/transfer Unit: 4.9 kg (11 lbs)

Mobile Platform: 5.4 kg (12 lbs)

Transfer Method

Centrifugal pump

Ordering information

Part #: 021-00000-00

Complete system includes storage Dewar, mobile platform, control/transfer unit, and 2.3-meter (7.5-ft) power cord



In keeping with a policy of ongoing product improvement, specifications are subject to change without notice.