

SPECIFICATIONS

Micromeritics® Subsieve AutoSizer (SAS) II



PHYSICAL

Height	55 cm (21.7 in)
Width	50 cm (19.7 in)
Depth	38 cm (15 in)
Weight	28 kg (61.7lbs)

ELECTRICAL

Voltage	120 - 240 VAC
Frequency	50 - 60 Hz
Current	1A

PARTICLE SIZE RANGE

Particle Size Range	0.2 - 75um
---------------------	------------

POROSITY RANGE

Porosity Range	0.2 - 0.9%
----------------	------------

COMPRESSION ACCURACY

Compression Accuracy	<0.05 mm
----------------------	----------

Due to continuous improvements, specifications are subject to change without notice.